

## DESCRIPTION

|                           |                                                                                                                                                                                                                                                                             |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Species Reactivity</b> | SARS-CoV-2                                                                                                                                                                                                                                                                  |
| <b>Specificity</b>        | Detects SARS-CoV-2 NSP8 in direct ELISAs.                                                                                                                                                                                                                                   |
| <b>Source</b>             | Monoclonal Mouse IgG <sub>2B</sub> Clone # 1044736                                                                                                                                                                                                                          |
| <b>Purification</b>       | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                  |
| <b>Immunogen</b>          | <i>E. coli</i> -derived SARS-CoV-2 NSP8 protein<br>Ala1-Gln198<br>Accession # YP_009725304.1                                                                                                                                                                                |
| <b>Conjugate</b>          | Alexa Fluor 350<br>Excitation Wavelength: 346 nm<br>Emission Wavelength: 442 nm                                                                                                                                                                                             |
| <b>Formulation</b>        | Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide<br><br>*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions. |

## APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

**Western Blot** Optimal dilution of this antibody should be experimentally determined.

## PREPARATION AND STORAGE

|                                |                                                                                                                   |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------|
| <b>Shipping</b>                | The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below. |
| <b>Stability &amp; Storage</b> | Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied                          |

## BACKGROUND

Non-structural protein 8 (NSP8) is one of several functional proteins released by ORF1a-encoded protease cleavage of the pp1a and pp1ab replicase polyproteins expressed from the coronavirus (CoV) genome (1). The NSPs are involved in the replication and transcription of the viral RNA and not incorporated within the virion particles. Coronaviruses include various highly pathogenic strains such as SARS-CoV, MERS-CoV and SARS-CoV2 that have had significant impact on humans as well as strains that have negatively impacted livestock. NSP8 is a small 198 amino acid protein that forms a dimer with NSP7 and subsequently assembles into a large hexadecameric structure (2). The NSP8 sequence is highly conserved across coronaviruses (2). The NSP8 monomers in each of two asymmetric units can adopt two different conformations: a "golf-club like" structure with the N-terminal shaft and C-terminal head or a bent N-terminal shaft and C-terminal head domain (2). The units stack to form a supercomplex like bricks with layers of NSP7 filling the spaces in between. The supercomplexes are stacked to form a channel with electrostatic properties that would allow RNA to pass through the channel, likely to facilitate efficient replication and transcription. NSP8 was also shown to be able to polymerize small oligomers in a sequence-specific fashion and was consequently proposed to act as an RNA primase for the viral RNA-dependent RNA polymerase (RdRp), NSP12, in SARS-CoV (3). In SARS-CoV-2, RdRp has been shown to have little activity without NSP8/7 acting as cofactors to form a complex (4) making NSP8 critical for viral polymerase activity. NSP8 was shown to interact with several other viral NSP proteins, including NSP2 and NSP9, (5) as well as multiple host cell proteins (6). NSP8 showed interaction with 3 different signal recognition particle (SRP) host cell proteins suggesting the virus hijacks the SEC61-mediated protein translocation pathway for entry into the endoplasmic reticulum and NSP8 may be involved in targeting (6).

## PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.