

## DESCRIPTION

|                                                                                                                                                                                                         |                                                                                             |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
| <b>Species Reactivity</b>                                                                                                                                                                               | SARS-CoV-2                                                                                  |
| <b>Specificity</b>                                                                                                                                                                                      | Detects SARS-CoV-2 NSP7 in direct ELISAs.                                                   |
| <b>Source</b>                                                                                                                                                                                           | Monoclonal Mouse IgG <sub>2A</sub> Clone # 1045511                                          |
| <b>Purification</b>                                                                                                                                                                                     | Protein A or G purified from hybridoma culture supernatant                                  |
| <b>Immunogen</b>                                                                                                                                                                                        | <i>E. coli</i> -derived SARS-CoV-2 NSP7 protein<br>Ser1-Gln83<br>Accession # YP_009725303.1 |
| <b>Conjugate</b>                                                                                                                                                                                        | Alexa Fluor 532<br>Excitation Wavelength: 534 nm<br>Emission Wavelength: 553 nm             |
| <b>Formulation</b>                                                                                                                                                                                      | Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide                                |
| *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions. |                                                                                             |

## APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

**Western Blot** Optimal dilution of this antibody should be experimentally determined.

## PREPARATION AND STORAGE

|                                |                                                                                                                   |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------|
| <b>Shipping</b>                | The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below. |
| <b>Stability &amp; Storage</b> | Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied                          |

## BACKGROUND

Non-structural protein 7 (NSP7) is one of several functional proteins released by ORF1a-encoded protease cleavage of the pp1a and pp1ab replicase polyproteins expressed from the coronavirus (CoV) genome (1). The NSPs are involved in the replication and transcription of the viral RNA and not incorporated within the virion particles. Coronaviruses include various highly pathogenic strains such as SARS-CoV, MERS-CoV and SARS-CoV-2 that have had significant impact on humans in addition to strains that have negatively impacted livestock. NSP7 is a small 83 amino acid protein that is highly conserved across coronaviruses (1). The NSP7 monomers are composed of a central core N-terminal tri-helical bundle with an additional short C-terminal helix (2). The monomeric units associate to form a large hexadecameric structure with NSP8 (2) where layers of NSP7 fill the spaces in between NSP8 units. The supercomplexes are stacked to form a channel with electrostatic properties that could allow RNA to pass through the channel to facilitate efficient replication and transcription. The NSP7/NSP8 supercomplex was thereby proposed to function as a primase for the viral RNA-dependent RNA polymerase (RdRp), NSP12 (3). In SARS-CoV-2, the RdRp has been shown to have little activity without NSP8/7 as cofactors (4) making NSP7/NSP8 critical for viral polymerase activity. NSP7 has also been shown to interact with several other viral NSP proteins, including NSP5, NSP9, and NSP13 (5) as well as multiple host cell proteins involved in membrane trafficking, signaling, and electron transport including potential drug targets such as COMT and PTGES (6).

## PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc., and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.