

## DESCRIPTION

|                           |                                                                                                                                                                                                                                                                             |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Species Reactivity</b> | Human                                                                                                                                                                                                                                                                       |
| <b>Specificity</b>        | Detects human CART aa 28-116 in direct ELISAs and Western blots. Also recognizes a truncated recombinant human CART containing amino acids 69-116.                                                                                                                          |
| <b>Source</b>             | Monoclonal Mouse IgG <sub>2B</sub> Clone # 113612                                                                                                                                                                                                                           |
| <b>Purification</b>       | Protein A or G purified from ascites                                                                                                                                                                                                                                        |
| <b>Immunogen</b>          | <i>E. coli</i> -derived recombinant human CART<br>Ala37-Leu116<br>Accession # Q16568                                                                                                                                                                                        |
| <b>Conjugate</b>          | Alexa Fluor 647<br>Excitation Wavelength: 650 nm<br>Emission Wavelength: 668 nm                                                                                                                                                                                             |
| <b>Formulation</b>        | Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide<br><br>*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions. |

## APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

|                             |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| <b>Western Blot</b>         | Optimal dilution of this antibody should be experimentally determined. |
| <b>Immunohistochemistry</b> | Optimal dilution of this antibody should be experimentally determined. |

## PREPARATION AND STORAGE

|                                |                                                                                                                   |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------|
| <b>Shipping</b>                | The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below. |
| <b>Stability &amp; Storage</b> | Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied                          |

## BACKGROUND

CART, also known as cocaine and amphetamine-related transcript, is a 10 kDa secreted polypeptide produced by neuronal cell types in a variety of locations. In human, the molecule is synthesized as a 116 amino acid (aa) precursor with a 27 aa signal sequence and an 89 aa mature segment. There are six C-terminal cysteines that form three intrachain disulfide bonds, and the potential exists for both sulfation and phosphorylation in the N-terminal region. There are numerous monobasic and dibasic sites for enzyme cleavage, and multiple, presumably proteolytically-generated short forms are known to exist. These include a 5K, 48 aa peptide (CART 42-89), and a 4.3K, 40 aa peptide (CART 49-89). Human CART is 98% aa identical to both mouse and rat CART (1-89). Although CART suppresses feeding and increases pain tolerance, there would appear to be differences between the various forms of CART on other functions. In addition, anatomical location is associated with differing lengths of CART. For example, 89 aa and 79 aa forms are found in the adrenal, while 48 aa and 44 aa forms are found in the hypothalamus.

## PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc., and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.