

## Human TrkA Alexa Fluor® 700-conjugated Antibody

Monoclonal Mouse IgG1 Clone # 165131 Catalog Number: FAB1751N

100 µg

DESCRIPTION	
Species Reactivity	Human
Specificity	Detects human TrkA in direct ELISAs. In direct ELISAs, this antibody does not cross-react with recombinant rat TrkA, recombinant human (rh) TrkB, or rhTrkC.
Source	Monoclonal Mouse IgG <sub>1</sub> Clone # 165131
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Mouse myeloma cell line NS0-derived recombinant human TrkA Ala33-Glu407 Accession # P04629
Conjugate	Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details.
	*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS	
Please Note: Optimal dilutions should be determined by each laboratory for each application	n. General Protocols are available in the Technical Information section on our website.
Recommended Concentration	Sample

Flow Cytometry	0.25-1 μg/10 <sup>6</sup> cells	See Below

DATA	
DATA Flow Cytometry	Detection of TrkA in Human Blood Monocytes by Flow Cytometry. Human peripheral blood monocytes were stained with Mouse Anti-Human CD14 PE-conjugated Monoclonal Antibody (Catalog # FAB3832P) and either (A) Mouse Anti-Human TrkA Alexa Fluor® 700-conjugated Monoclonal Antibody (Catalog # FAB1751N) or (B) Mouse IgG <sub>1</sub> Alexa Fluor 700 Isotype Control (Catalog # IC002N). View our protocol for Staining Membrane- associated Proteins.
(D) AND S	TOPACE
PREPARATION AND S Shipping	TORAGE The product is shipped with polar packs. Upon receipt, st
Stability & Storage	Protect from light. Do not freeze.
- •	<ul> <li>12 months from date of receipt 2 to 8 °C as supp</li> </ul>

### BACKGROUND

TrkA, the product of the proto-oncogene trk, is a member of the neurotrophic tyrosine kinase receptor family that has three members. TrkA, TrkB, and TrkC. TrkA binds NGF and NT-3, while TrKB binds NT-4 and BDNF, and TrkC binds NT-3. All Trk family proteins share a conserved complex subdomain organization consisting of two cysteine-rich segments, a cluster of three leucine-rich motifs, and two immunoglobulin-like domains in the extracellular region, coupled to an intracellular region that contains a tyrosine kinase domain. Two distinct TrkA isoforms that differ by virtue of a 6-amino acid (aa) insertion after Pro392 in their extracellular domain have been identified. The longer TrkA isoform is the only isoform expressed within neuronal tissues, whereas the shorter TrkA is expressed mainly in non-neuronal tissues. A third TrkA variant is expressed by neural stem cells. This isoform shows a Val substitution for aa 194-283 and an absence of the aforementioned 6 aa insert. NGF binds to TrkA with low affinity and activates its cytoplasmic kinase, initiating a signaling cascade that mediates neuronal survival and differentiation. Higher affinity binding of NGF requires the coexpression of TrkA with the p75 NGF receptor (NGF R), a member of the tumor necrosis factor receptor superfamily. NGF R binds all neurotrophins with low affinity and modulates Trk activity as well as alters the specificity of Trk receptors for their ligands. NGF R can also mediate cell death when expressed independent of Trk. Over aa 33-407, human TrkA shares 79% aa sequence identity with mouse TrkA.

## Rev. 5/20/2019 Page 1 of 2

biotechne

Global bio-techne.com info@bio-techne.com techsupport@bio-techne.com TEL +1 612 379 2956 USA TEL 800 343 7475 Canada TEL 855 668 8722 China TEL +86 (21) 52380373 Europe | Middle East | Africa TEL +44 (0)1235 529449



# Human TrkA Alexa Fluor® 700-conjugated Antibody

Monoclonal Mouse IgG<sub>1</sub> Clone # 165131 Catalog Number: FAB1751N 100 μg

### PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Rev. 5/20/2019 Page 2 of 2



**Global** bio-techne.com info@bio-techne.com techsupport@bio-techne.com TEL +1 612 379 2956 USA TEL 800 343 7475 Canada TEL 855 668 8722 China TEL +86 (21) 52380373 Europe | Middle East | Africa TEL +44 (0)1235 529449