

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human LILRB3/CD85a/ILT5 in direct ELISA.
Source	Monoclonal Mouse IgG ₁ Clone # 1057017
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Human embryonic kidney cell, HEK293-derived human LILRB3/CD85a/ILT5 Gly24-Glu443 Accession # O75022
Conjugate	Alexa Fluor 350 Excitation Wavelength: 346 nm Emission Wavelength: 442 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Flow Cytometry	Titration recommended for optimal concentration with starting range of 0.1-1 µg/1 million cells. Sample used for this experiment was PBMCs with CD14 costain.
-----------------------	---

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

Leukocyte immunoglobulin-like receptor subfamily B (LILRB3), also known as ILT5, LIR3, and CD85a, is an immunoglobulin superfamily member that is involved in immune regulation. Subfamily B members have cytoplasmic immunoreceptor tyrosine-based inhibitory motifs (ITIMs) that inhibit signaling events via phosphatase SHP-1. Subfamily A members are activating receptors that lack ITIMs and signal through association with Fcγ₁ (1, 2). Mature LILRB3 is a highly polymorphic 85-95 kDa glycoprotein that consists of a 420 amino acid (aa) extracellular domain (ECD) with four Ig-like domains, a 21 aa transmembrane segment, and a 167 aa cytoplasmic domain with three ITIMs (3). Alternative splicing generates an isoform with a 17 aa insertion in the juxtamembrane ECD. In mouse and rat, the LILRB3 gene encodes the PIR-B protein which has six Ig-like domains. Rodent PIR-B and human LILRB3 share 55% aa sequence identity within common regions of their ECDs. LILRB3 is expressed on the surface of peripheral monocytes, neutrophils, eosinophils, basophils, and mast cell progenitors (4-6). Triggering of LILRB3 inhibits the activation of macrophages, mast cells, neutrophils, basophils, and B cells (5, 7). On osteoclast precursors, LILRB3 ligation inhibits RANK L/TRANCE or M-CSF induced differentiation (8). LILRB3 can also bind to ligands exposed on necrotic tumor cells (9). Both PIR-B and LILRB3 are receptors for S. aureus, and activation of these receptors by bacteria influences the innate immune response triggered by TLRs (3). R&D Systems in-house testing indicates that LILRB3 binds to Angiopoietin-like 7, consistent with the demonstrated functional interactions between other members of these protein families (10). In the mouse CNS, PIR-B functions as a receptor for the myelin proteins Nogo, MAG, and OMgp and mediates their inhibitory action on neurite outgrowth and axon regeneration (11). Upon binding to MAG, PIR-B associates with TrkB and NGF R/p75 in cerebellar granule neurons (12).

References:

1. Thomas, R. *et al.* (2010) Clin. Rev. Allergy Immunol. **38**:159.
2. Anderson, K.J. and R.L. Allen (2009) Immunology **127**:8.
3. Nakayama, M. *et al.* (2007) J. Immunol. **178**:4250.
4. Tedla, N. *et al.* (2003) Proc. Natl. Acad. Sci. USA **100**:1174.
5. Sloane, D.E. *et al.* (2004) Blood **104**:2832.
6. Tedla, N. *et al.* (2008) J. Leukoc. Biol. **83**:334.
7. Uehara, T. *et al.* (2001) J. Clin. Invest. **108**:1041.
8. Mori, Y. *et al.* (2008) J. Immunol. **181**:4742.
9. Jones, D.C. *et al.* (2016) Oncotarget PMID 26769854.
10. Zheng, J. *et al.* (2012) Nature **485**:656.
11. Atwal, J.K. *et al.* (2008) Science **322**:967.
12. Fujita, Y. *et al.* (2011) Cell Death Dis. **2**:e198.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.