

DESCRIPTION

Species Reactivity	Human
Specificity	Detects recombinant human BLAME in ELISAs.
Source	Monoclonal Mouse IgG _{2A} Clone # 250020
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Mouse myeloma cell line NS0-derived recombinant human BLAME/SLAMF8 Ala23-Asp233 Accession # Q9P0V8
Conjugate	Alexa Fluor 647 Excitation Wavelength: 650 nm Emission Wavelength: 668 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

ELISA Capture (Matched Antibody Pair)	Optimal dilution of this antibody should be experimentally determined.
ELISA Detection (Matched Antibody Pair)	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

BLAME (B-lymphocyte activator macrophage expressed), also known as SLAM family member 8, is a type I transmembrane protein that belongs to the CD2 subset of immunoglobulin superfamily cell receptors. CD2 family proteins function as adhesion molecules and modulators of immune responses (1, 2). Mature human BLAME consists of a 211 amino acid (aa) ECD that contains two Ig V-like domains, a 21 aa transmembrane segment, and a 31 aa cytoplasmic tail that lacks recognizable signaling motifs (3). Within the ECD, human BLAME shares 19% - 26% aa sequence identity with human 2B4, CD2, CD2F-10, CD48, CD58, CD84, CD229, CRACC, NTB-A, and SLAM. It shares 79% aa sequence identity with the ECD of mouse BLAME. BLAME is expressed on dendritic cells and IFN-γ stimulated monocytes. Overexpression of BLAME in bone marrow cells leads to an increase in the peritoneal B1b population of B lymphocytes (3).

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.