

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human CRIM1 in direct ELISAs and Western blots.
Source	Monoclonal Mouse IgG ₁ Clone # 253511
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Mouse myeloma cell line NS0-derived recombinant human CRIM1 Leu35-Ser939 Accession # Q9NZV1
Conjugate	Alexa Fluor 488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

Cysteine rich motor neuron 1 (CRIM1) is a type I transmembrane glycoprotein of the chordin-like cysteine-rich repeat (CRR) family of BMP inhibitors (1 - 4). The ~130 kDa, 1036 amino acid (aa) CRIM1 contains a 34 aa signal sequence, a 905 aa extracellular domain (ECD), a 21 aa transmembrane domain and a 76 aa cytoplasmic domain. The ECD includes an N-terminal IGF-binding protein-like motif and six chordin-like von Willebrand C-type CRRs. The ECD can be released from the cell, presumably by proteolytic processing (4). Human CRIM1 ECD shows 88%, 88%, 91%, 86%, 87%, 83%, and 72% aa identity with mouse, rat, dog, cow, opossum, chick, and zebrafish CRIM1 ECD, respectively. CRIM1 can interact with TGF-β family ligands, including BMPs 2, 4, and 7, via its CRR domains (4). It binds BMPs intracellularly and antagonizes them by lowering their expression, processing and secretion (4). CRIM1 is expressed in the developing spinal cord in the floor plate and developing motor neurons (1). It is also expressed by perivascular smooth muscle cells and aligns at points of cell-cell contact during endothelial cell capillary formation (2). Endothelial cell expression in vitro appears to be specific to cells that are adherent and growing (2). CRIM1 is also expressed in a spatially and temporally restricted manner in the developing lens, limbs, kidney, teeth, and testis (5). Studies where CRIM1 expression is manipulated in developing mouse, chick, and zebrafish support its involvement in regulation of vascular and somitic development and organogenesis (5 - 7).

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.