

DESCRIPTION


Species Reactivity	Human
Specificity	Detects human DEC-205/CD205 in direct ELISAs. In direct ELISAs, no cross-reactivity with recombinant mouse DEC-205 is observed.
Source	Monoclonal Mouse IgG ₁ Clone # 523203
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant human DEC-205/CD205 Cys216-Cys501 Accession # O60449
Conjugate	Fluorescein Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm (FITC)
Formulation	Supplied in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Flow Cytometry	10 μ L/10 ⁶ cells	See Below

DATA



PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. <ul style="list-style-type: none"> ● 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

DEC-205, also known as CD205 and lymphocyte antigen 75 (Ly 75), is a type I transmembrane protein that is primarily expressed on dendritic cells and thymic epithelial cells. The extracellular region of DEC-205 contains ten C-type lectin-like domains, a fibronectin type II domain and a ricin B-type lectin domain. DEC-205 functions as an endocytic receptor for antigens. The recombinant protein used to generate the anti-human DEC-205 antibody contains the first two C-type lectin domains.