

DESCRIPTION

Species Reactivity	Mouse
Specificity	Detects mouse Neuropoietin/NP in direct ELISAs and Western blots.
Source	Monoclonal Rat IgG ₁ Clone # 321525
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant mouse Neuropoietin/NP Ala23-Ala204 Accession # P83714
Conjugate	Alexa Fluor 488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

Neuropoietin (NP; also known as cardiotrophin-2) is a 22 kDa member of the IL-6 family of cytokines. Considered to be the product of a gene duplication event involving cardiotrophin-1 (CT-1), it helps to define a subfamily within the IL-6 family that includes CT-1, CLC and CTNf. Mouse neuropoietin is synthesized as a 204 amino acid (aa) precursor that contains a 22 aa signal sequence and a 192 aa mature segment. The secreted molecule is characterized by the presence of four α-helices, typical of hematopoietic superfamily molecules. Mature mouse neuropoietin shares 88%, 90% and 95% aa identity to chimpanzee, canine and rat neuropoietin, respectively. The human gene is suggested to have evolved towards a pseudogene, a point of interest in that neuropoietin is reported to signal through the CNTF complex (*i.e.*, gp130, CNTF Rα and LIF R). NP will mediate motor neuron survival, and appears to be selectively expressed in the embryo by tissues involved with nervous system development (1).

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.