

Mouse Integrin β 4/CD104 Alexa Fluor® 700-conjugated Antibody

Monoclonal Rat IgG_{2A} Clone # 308601

Catalog Number: FAB4054N

100 μ g

DESCRIPTION

Species Reactivity	Mouse
Specificity	Detects mouse Integrin β 4/CD104 in direct ELISAs and Western blots. Does not cross-react with recombinant mouse Integrin β 1, β 2, β 6, β 7, or recombinant human Integrin β 3.
Source	Monoclonal Rat IgG _{2A} Clone # 308601
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Chinese hamster ovary cell line CHO-derived recombinant mouse Integrin β 4/CD104 Asn29-Ser711 Accession # NP_001005608
Conjugate	Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

	Recommended Concentration	Sample
Flow Cytometry	0.25-1 μ g/10 ⁶ cells	CMT-93 mouse rectal carcinoma cell line

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

Integrin β 4, also known as CD104, is a 200 kDa type I transmembrane (TM) glycoprotein that associates only with Integrin α 6/CD29. Integrin α 6 β 4 is predominantly expressed by epithelial cells, binds laminins, and is essential for formation of hemidesmosomes and connection of the dermis to the epidermis. Mouse integrin β 4 contains a 680 amino acid (aa) extracellular domain (ECD) with a metal-binding site and four cysteine-rich repeats, a 22 aa TM segment and an unusually long 1073 aa cytoplasmic tail that organizes hemidesmosome components. Isoforms with shortened cytoplasmic tails have been described. Mouse and human β 4 ECD show 88% aa identity.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.