

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human ADAMTS4 in direct ELISAs and Western blots. In direct ELISAs and Western blots, no cross-reactivity with recombinant human (rh) ADAMTS1, rhADAMTS5, rhADAMTS13, or rhADAMTSL1.2 is observed.
Source	Monoclonal Mouse IgG _{2A} Clone # 416610
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Mouse myeloma cell line NS0-derived recombinant human ADAMTS4 Phe213-Cys685 Accession # O75173
Conjugate	Alexa Fluor 488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

Neutralization	Optimal dilution of this antibody should be experimentally determined.
Western Blot	Optimal dilution of this antibody should be experimentally determined.
Immunocytochemistry	Optimal dilution of this antibody should be experimentally determined.
Immunoprecipitation	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

ADAMTS4 (a disintegrin and metalloproteinase with thrombospondin motifs 4), also known as aggrecanase-1, is a member of the family of secreted zinc proteases with a multi-domain structure (1-3). The protein precursors consist of a signal peptide and the following domains: pro, catalytic, disintegrin-like, TS type 1 motif, cysteine-rich, and spacer. It is the only ADAMTS identified that has one TS type I motif. It is an active protease effectively cleaving α-2-macroglobulin and aggrecan at multiple sites, and is inhibited by TIMP-3 with inhibition constants in subnanomolar range (4-6). It receives great attention due to the elevation in its mRNA level after treatment with Interleukin-1 (7). However, in a mouse model of osteoarthritis, ADAMTS4 knock-out mice did not exhibit any significant protective effect (8). ADAMTS4 consists of a signal peptide (aa 1 to 51), a pro domain (aa 52 to 212) and a mature chain (aa 213 to 837) containing the following domains: catalytic (aa 218 to 425), disintegrin (aa 437 to 519), TSP type-1 (aa 520 to 575), Cys-rich (aa 577 to 685) and spacer (aa 686 to 837). The amino acid sequence of rhADAMTS4 is 100%, 97% and 94% identical to that of chimpanzee, dog, and mouse/rat/bovine.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc., and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.