

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human S100A7 in direct ELISAs. In direct ELISAs, no cross-reactivity with recombinant human (rh) S100A1, A2, A4, A6, A8, A9, A10, A11, A14, A16, rhS100B, rhS100P, recombinant mouse S100A1, A4, A6, A8, A9, A10, A11, A13, or A16 is observed.
Source	Monoclonal Mouse IgG ₁ Clone # 577513
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant human S100A7 Met1-Gln101 Accession # P31151
Conjugate	Alexa Fluor 488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

Immunohistochemistry Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

S100A7 (also Psoriasin) is an 11-13 kDa member of the S101 family, EF-hand superfamily of Ca-binding proteins. It is produced by keratinocytes and exists as both an intracellular and extracellular polypeptide. Intercellularly, it binds to E/FABP and may participate in fatty acid metabolism. Extracellularly, it is both monomeric and heterooligomeric and exhibits antibacterial properties. Human S100A7 is 101 amino acids (aa) in length. It contains two EF-hand motifs (aa 13-48 and 50-85), three zinc-binding His residues, one high-affinity Ca-binding site (aa 63-74) and an overlapping antibacterial sequence (aa 35-80). There appear to be no meaningful rodent structural orthologs. Human S100A7 is 93% aa identical to human S100A15/7A, and the S100A7 gene has a duplication (S100A7B) that is 54% aa identical to S100A7.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.