

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human Desmocollin-1 in direct ELISA.
Source	Monoclonal Rabbit IgG Clone # 2906A
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Mouse myeloma cell line, NS0-derived human Desmocollin-1 Arg135-Asn686 Accession # Q08554
Conjugate	Alexa Fluor 594 Excitation Wavelength: 590 nm Emission Wavelength: 617 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Immunohistochemistry Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

Desmocollin-1 (DSC1) is a 100 kDa-110 kDa transmembrane glycoprotein in the cadherin family of calcium dependent adhesion molecules (1-3). Human DSC1 is synthesized with a 105 amino acid (aa) propeptide. The mature protein consists of a 557 aa extracellular domain (ECD) that contains five cadherin-like domains, a 23 aa transmembrane segment, and a 180 aa cytoplasmic domain (4, 5). Within the ECD, human DSC1 shares 79% and 82% aa sequence identity with bovine and mouse DSC1, respectively, and 53% with human DSC2 and DSC3. An alternately spliced isoform has a substituted and truncated cytoplasmic domain (5). DSC1 is one of the principal components of desmosomes which form adhesive contacts between epithelial cells (1, 2). It is predominantly expressed in the outer layers of stratified epithelia of the skin, tongue, and hair follicle root sheath (6, 7). It is required for both promoting epidermal differentiation and maintaining epidermal strength (8). DSC2 and DSC3, by contrast, are preferentially localized in the basal and suprabasal layers of the epidermis (1). DSC1 is not expressed in normal colon epithelium, but it is induced during colon carcinogenesis (9). Loss of DSC1 function in the skin disorders IgA pemphigus and Netherton Syndrome results from autoimmune targeting or enhanced KLK-7 mediated proteolysis, respectively (10-12).

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.