

DESCRIPTION

Species Reactivity	Human/Rat
Specificity	Detects human and rat DYRK2 in Western blots.
Source	Monoclonal Mouse IgG _{2B} Clone # 599542
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant human DYRK2 aa 1-107 Accession # NP_003574
Conjugate	Alexa Fluor 350 Excitation Wavelength: 346 nm Emission Wavelength: 442 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot	Optimal dilution of this antibody should be experimentally determined.
Immunohistochemistry	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

DYRK2 (Dual-specificity tyrosine [Y] phosphorylation regulated kinase 2) is a 58-62 kDa member of the MNB/DYRK subfamily, CMGC Ser/Thr protein kinase family of enzymes. It is expressed in testis and shows dual substrate specificity; autophosphorylation on Tyr382 to self-activate, and a Ser/Thr phosphorylation of target molecules. Substrates include NFATc, glycogen synthase, and p53. p53 phosphorylation on Ser46 initiates cell apoptosis. Human DYRK2 is 601 amino acids (aa) in length. It contains one kinase catalytic domain (aa 222-535). One isoform variant shows an alternate start site at Met74. This shorter isoform was used for immunization. Thus, over aa 1-107 (or aa 74-180 of the long form), human DYRK2 shares 93% aa identity with mouse DYRK2.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.