

DESCRIPTION

Species Reactivity	Rat
Specificity	Detects rat Agrin in direct ELISAs and Western blots.
Source	Monoclonal Mouse IgG ₁ Clone # 76805
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>S. frugiperda</i> insect ovarian cell line Sf21-derived recombinant rat Agrin Ala1153-Pro1959 (Pro1788-Ser1798 del) Accession # P25304
Conjugate	Alexa Fluor 750 Excitation Wavelength: 749 nm Emission Wavelength: 775 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

Agrin is a 400-600 kDa heparan sulfate proteoglycan component of the extracellular matrix. The N-terminal half of rat Agrin, which mediates ECM interactions, contains nine Kazal-type protease inhibitor domains, two Laminin EGF-like domains, and one SEA domain. The C-terminal half contains four EGF-like repeats and three Laminin globular G domains. Human Agrin also contains a Laminin-binding N-terminal Agrin domain (NtA), and mouse and chick Agrin include the NtA domain only by the use of an alternate promoter. Additional isoforms are generated by alternate splicing at sites Y and Z in the C-terminal half of rat Agrin (known as A and B, respectively in chick). Agrin isoforms that contain an insert at site Z (Z+ forms) are known as neural Agrin and are selectively produced by motoneurons. Other isoforms are known as muscle Agrin and are additionally expressed in non-neuronal tissues, particularly in basement membranes of the lung and kidney (1-3). This recombinant protein consists of the C-terminal half of rat Agrin and contains a nine amino acid (aa) insert at the Z site. It shares 59%, 80%, and 94% aa sequence identity with comparable regions of chick, human, and mouse Agrin, respectively. The C-terminal half of Z- and Z+ Agrin binds to α-Dystroglycan and mediates adhesion between motoneurons and myotubes at the neuromuscular junction (NMJ) (4-6). In contrast, only Z+ Agrin is effective at inducing clustering of the postsynaptic Acetylcholine Receptor (AChR) and presynaptic motoneuron differentiation (7, 8). Agrin-induced AChR clustering requires a myotube receptor complex that contains α-Dystroglycan, MuSK, and LRP4 (4, 9-11). Agrin exhibits many functions in addition to NMJ development. It is enriched in senile Alzheimer's disease plaques where it binds the Aβ (1-40) peptide and promotes amyloid fibril formation (12). It regulates neuronal excitability by binding and inhibiting the α3 subunit of the neuronal Na/K ATPase (13). It functions as an epithelial cell attachment receptor for HIV-1 through interactions with the gp41 coat protein (14). During T cell activation, Agrin contributes to formation of the immunological synapse and regulates the threshold of T cell activation (15).

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.