

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human ELF3 in direct ELISAs. In direct ELISAs, no cross-reactivity with recombinant human ELF5 is observed.
Source	Monoclonal Mouse IgG ₁ Clone # 662516
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant human ELF3 Met1-Gly173 Accession # P78245
Conjugate	Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Immunohistochemistry Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

ELF3 (E74-Like Factor 3; also ESE-1, ESX, ERT and JEN) is a 41-43 kDa member of the ETS family of proteins. During uncomplicated (non-inflammatory) periods of cell differentiation, ELF3 is expressed exclusively by epithelial cells, repressing genes needed during early differentiation, and promoting genes needed for full differentiation. Under conditions of inflammation, cells such as monocytes, endothelial cells and chondrocytes will express ELF3 and produce molecules such as Ang1 and COX2. Human ELF3 is 371 amino acids (aa) in length. It contains one PNT/pointed dimerization domain (aa 46-132), a protein stabilizing PEST sequence (aa 210-225), an A/T Hook region that binds to AT-rich DNA sequences (aa 236-252), and an ETS DNA binding domain (aa 273-355). ELF3 interacts with CREBBP, EP300, KU70 and KU86. There is one splice variant that shows a deletion of aa 174-200. Over aa 1-173, human ELF3 shares 87% aa identity with mouse ELF3.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.