

DESCRIPTION

Species Reactivity	Human
Specificity	Detects Fatty Acid Synthase/FASN in direct ELISAs and Western blots. In direct ELISAs, approximately 25% cross-reactivity with recombinant mouse (rm) Ephrin-A2 is observed and no cross-reactivity with recombinant human (rh) Ephrin
Source	Monoclonal Mouse IgG _{2B} Clone # 647114
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant human Fatty Acid Synthase/FASN Ser10-Cys212 Accession # P49327
Conjugate	Alexa Fluor 647 Excitation Wavelength: 650 nm Emission Wavelength: 668 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

FASN (Fatty Acid Synthase; also FAS) is a cytosolic 270-280 kDa member of the short chain dehydrogenase/reductase family 27X. It is expressed in neurons, skeletal muscle, adipocytes and hepatocytes, and both catalyzes the formation of palmitic acid from acetylCoA and malonylCoA, and likely mediates the transfer of fatty acids to peptides. Human FASN is 2511 amino acids (aa) in length and contains multiple domains, including a β -ketoacyl synthase domain (aa 1-414), an acyl and malonyl transferase domain (aa 429-817), an enoyl reductase domain (aa 1563-1863), a β -ketoacyl reductase domain (aa 1864-2118) and a thioesterase domain (aa 2207-2511). FASN exists as an antiparallel homodimer. There is one theoretical alternative start site at Met2073. Over aa 9-212, human FASN shares 92% aa identity with mouse FASN.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.