

Human TSHR Alexa Fluor® 700-conjugated Antibody

Monoclonal Mouse IgG₁ Clone # 484404
Catalog Number: FAB65342N
100 µg

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human TSHR in direct ELISAs.
Source	Monoclonal Mouse IgG ₁ Clone # 484404
Purification	Protein A or G purified from ascites
Immunogen	NS0 mouse myeloma cell line transfected with human TSHR Accession # P16473
Conjugate	Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the Technical Information section on our website.

	Recommended Concentration	Sample
Flow Cytometry	0.25-1 µg/10 ⁶ cells	HEK293 Human Cell Line Transfected with Human TSH R and eGFP

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. • 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

Thyroid stimulating hormone receptor (TSH R; also thyrotropin receptor) is an 85 kDa (unglycosylated) member of the G-protein coupled receptor 1 family. Human TSH R is synthesized as a 764 amino acid (aa) precursor that contains a 20 aa signal sequence, a 393 aa extracellular domain (ECD), a 269 aa membrane spanning domain consisting of seven transmembrane segments, and an 82 aa cytoplasmic domain. The ECD contains six leucine-rich repeats (LRR) and six potential sites of N-linked glycosylation. There are two isoforms, long and short, produced by a splicing variant corresponding to aa 253-764 in the long isoform that is missing in the short isoform. Human TSH R is 86-87% aa identical to mouse and rat TSH R. TSH R is the receptor for thyroid stimulating hormone, which plays a central role in controlling thyroid cell metabolism.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc., and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.