

Mouse Semaphorin 4F Alexa Fluor® 750-conjugated Antibody

Monoclonal Rat IgG_{2A} Clone # 780225

Catalog Number: FAB7200S

100 µg

DESCRIPTION

Species Reactivity	Mouse
Specificity	Detects mouse Sema-4F in ELISAs.
Source	Monoclonal Rat IgG _{2A} Clone # 780225
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Mouse myeloma cell line NS0-derived recombinant mouse Sema-4F Arg41-Gly667 Accession # Q9Z123
Conjugate	Alexa Fluor 750 Excitation Wavelength: 749 nm Emission Wavelength: 775 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Flow Cytometry	0.25-1 µg/10 ⁶ cells	Mouse cortical stem cells and bEnd.3 mouse endothelioma cell line

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. <ul style="list-style-type: none"> 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

Semaphorin 4F (Sema4F; also Sema W) is a 90-100 kDa, class IV member of the semaphorin family of proteins. It is expressed by select cell types, including Schwann cells and neurons, and likely acts as an inhibitor of cell proliferation by Schwann cells, and a regulator of postsynaptic glutamatergic synapses on neuron dendrites. Mature mouse Sema4F is a type I transmembrane protein that is 737 amino acids (aa) in length. It contains a 627 aa extracellular region (aa 41-667) that is characterized by the presence of one Sema domain (aa 71-493), a PSI region (aa 518-569), and a C2-type Ig-like domain (aa 586-641). Over aa 41-667, mouse Sema4F shares 92% and 98% aa identity with human and rat Sema4F, respectively.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.