

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human ZEB2/SIP1 in direct ELISAs.
Source	Monoclonal Mouse IgG _{2B} Clone # 923328
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant human ZEB2/SIP1 Asn363-Lys537 Accession # O60315
Conjugate	Alexa Fluor 750 Excitation Wavelength: 749 nm Emission Wavelength: 775 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Intracellular Staining by Flow Cytometry	0.25-1 µg/10 ⁶ cells	A549 human lung carcinoma cell line fixed and permeabilized with FlowX FoxP3 Fixation & Permeabilization Buffer Kit (Catalog # FC012)

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. <ul style="list-style-type: none"> 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

SIP1 (Smad-interacting protein 1), also called ZEB2 (Zinc finger E-box-binding homeobox 2), is a nuclear transcription factor. SIP1 contains seven C2H2-type zinc finger domains within amino acids (aa) 211-334 and 999-1076, a Smad-MHZ binding domain (aa 437-487), an atypical homeobox domain (aa 644-703), and phosphorylation, sumoylation and acetylation sites. The 1214 aa SIP1 gives a calculated molecular weight of 136 kDa, but may actually appear closer to 200 kDa due to modifications. A 1190 aa isoform lacks aa 111-134. Within aa 363-537, human SIP1 shares 98% and 97% aa sequence identity with mouse and rat SIP1, respectively. SIP1 is highly expressed in postmitotic neocortical cells and influences cell fate decisions in embryonic brain development. Point mutations causing underexpression of SIP1 are associated with Mowat-Wilson syndrome (MWIS), also known as Hirschprung disease mental retardation syndrome.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.