

Mouse IL-38/IL-1F10 Alexa Fluor® 350-conjugated Antibody

Monoclonal Rat IgG_{2A} Clone # 798036R

Catalog Number: FAB77741U

100 µg

DESCRIPTION

Species Reactivity	Mouse
Specificity	Detects mouse IL-38/IL-1F10 in direct ELISAs and Western blots.
Source	Monoclonal Rat IgG _{2A} Clone # 798036R
Purification	Protein A or G purified from cell culture supernatant
Immunogen	<i>E.coli</i> -derived recombinant mouse IL-38/IL-1F10 Met1-Arg152 Accession # Q8R459
Conjugate	Alexa Fluor 350 Excitation Wavelength: 346 nm Emission Wavelength: 442 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

	Recommended Concentration	Sample
Flow Cytometry	0.25-1 µg/10 ⁶ cells	RAW 264.7 mouse monocyte/macrophage cell line fixed with Flow Cytometry Fixation Buffer (Catalog # FC004) and permeabilized with Flow Cytometry Permeabilization/Wash Buffer I (Catalog # FC005)

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

IL-1F10, also known as IL-1HY2, is a secreted protein that shares sequence and structural similarity with IL-1 receptor antagonist and other IL-1 family members. It is expressed in basal epithelia of the skin and in proliferating germinal center CD20⁺ B cells. IL-1F10 binds soluble IL-1RI in vitro and may participate in regulation of immune responses. Mouse IL-1F10 shares 82% and 93% aa identity with human and rat IL-1F10, respectively. IL-1F10 is predicted to be 17 kDa.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc., and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.