

## DESCRIPTION

<b>Species Reactivity</b>	Mouse
<b>Specificity</b>	Detects mouse TGN38 in ELISAs.
<b>Source</b>	Monoclonal Rat IgG <sub>2B</sub> Clone # 831629
<b>Purification</b>	Protein A or G purified from hybridoma culture supernatant
<b>Immunogen</b>	<i>E. coli</i> -derived recombinant mouse TGN38 His320-Leu349 Accession # Q62313
<b>Conjugate</b>	Alexa Fluor 405 Excitation Wavelength: 405 nm Emission Wavelength: 421 nm
<b>Formulation</b>	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

## APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

<b>Western Blot</b>	Optimal dilution of this antibody should be experimentally determined.
<b>Immunocytochemistry</b>	Optimal dilution of this antibody should be experimentally determined.

## PREPARATION AND STORAGE

<b>Shipping</b>	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
<b>Stability &amp; Storage</b>	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

## BACKGROUND

TGN38A (also Tgoln1, Trans-Golgi network integral membrane protein 1) is an integral membrane protein associated with intracellular protein trafficking. Although its predicted MW is 38 kDa, due to a high proline content and extensive polysialylation, it runs anomalously at 80-90 kDa in SDS-PAGE. Tgoln1 is ubiquitously expressed and is generally associated with the trans-Golgi complex, a structure adjacent to the trans-component of the Golgi Apparatus. In this complex, both secretory and membrane proteins are sorted and forwarded to various compartments such as lysosomes, endosomes, the cell membrane and secretory granules. TGN38A apparently cycles between the trans-Golgi network and the cell membrane, returning to the trans-Golgi via the endosomal system. Part of its mobility is mediated by an interaction between TGN38A, neurabin and actin. Integrin β1 is one molecule suggested to be transported by TGN38A. Mature mouse TGN38A is a 336 amino acid (aa) type I transmembrane glycoprotein. It contains a 281 aa extracellular region (aa 18-298) plus a 34 aa cytoplasmic domain. The extracellular region possesses six sequential octapeptide repeats (aa 131-178) plus two utilized phosphorylation sites (Ser230 and Ser231). The cytoplasmic domain contains an endocytosis signal (Ser344-Leu349) that requires a free hydroxyl on Ser344 for proper routing. The mouse genome has two Tgn38 genes, Tgn38A and Tgn38B/Tgoln2. They differ in two ways. First there is a two aa insertion after Pro47, and an eight aa insertion after Thr154 in Tgn38B and, second, Tgn38B expression is restricted, occurring in mouse strain ICR, while Tgn38A is widespread and found in strains ICR, BALB/c and C57BL/6. Over aa 320-349, mouse and rat Tgoln1 are identical in aa sequence. There is no strict human structural counterpart to mouse Tgoln1/Tgn38A. Human does, however, possess a similar functioning molecule termed TGN46/48/51, and over the short aa stretch 320-349, mouse and human TGN molecules share 96% aa sequence identity.

## PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.