

DESCRIPTION

Species Reactivity	Human/Mouse/Rat
Specificity	Detects human Histone H2AY in direct ELISAs and Western blots.
Source	Monoclonal Mouse IgG ₁ Clone # 866416
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant human Histone H2AY Gly230-Thr348 Accession # O75367
Conjugate	Alexa Fluor 750 Excitation Wavelength: 749 nm Emission Wavelength: 775 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot	Optimal dilution of this antibody should be experimentally determined.
Immunocytochemistry	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

Histone H2AY, also known as macro H2A1 (MH2A1) and H2AFY, is a 40 kDa variant of the histone family, basic nuclear proteins responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Ubiquitously expressed, H2AY replaces conventional H2A histones in a subset of nucleosomes, where it represses transcription and participates in stable X chromosome inactivation. H2AY expression is also enriched in senescence-associated heterochromatin. Alternative splicing creates multiple transcript variants, encoding three distinct isoforms. H2AY expressed in blood has been identified as a potential biomarker associated with the activity and therapeutic response of Huntington disease, an autosomal dominant neurodegenerative disorder. Over amino acids 230-348, human H2AY shares 98% identity with mouse and rat H2AY.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.