

## DESCRIPTION


<b>Species Reactivity</b>	Mouse
<b>Specificity</b>	Detects mouse DC-SIGN/CD209 in direct ELISAs.
<b>Source</b>	Recombinant Monoclonal Rat IgG <sub>2B</sub> Clone # 902404
<b>Purification</b>	Protein A or G purified from cell culture supernatant
<b>Immunogen</b>	NS0 mouse myeloma cell line transfected with mouse DC-SIGN/CD209 Met1-Lys238 Accession # Q91ZX1
<b>Conjugate</b>	Phycoerythrin Excitation Wavelength: 488 nm Emission Wavelength: 565-605 nm
<b>Formulation</b>	Supplied in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details.  *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

## APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
<b>Flow Cytometry</b>	10 $\mu$ L/10 <sup>6</sup> cells	See Below

## DATA



## PREPARATION AND STORAGE

<b>Shipping</b>	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
<b>Stability &amp; Storage</b>	<b>Protect from light. Do not freeze.</b> <ul style="list-style-type: none"> <li>12 months from date of receipt, 2 to 8 °C as supplied.</li> </ul>

## BACKGROUND

Mouse Dendritic Cell-specific ICAM-3 Grabbing Non-integrin (DC-SIGN)/CD209, also known as CD209 Antigen-like Protein A, is a member of the C-type lectin family (1). Mouse DC-SIGN/CD209 is a 33 kDa, 238 amino acid (aa) type II transmembrane protein (2). The extracellular region contains a Ca<sup>2+</sup>-dependent carbohydrate-binding lectin domain (2). In addition to the full-length mouse DC-SIGN/CD209, a second, truncated splice variant has been reported. DC-SIGN/CD209 is not well conserved between mouse and human with the extracellular domain sharing only 63% aa identity. The DC-SIGN/CD209 lectin domain binds mannose oligosaccharides on pathogens including HIV as well as self glycoproteins including ICAMs (2-4). DC-SIGN/CD209 is expressed on dendritic cells (DC) and inflammatory macrophages and contributes to antigen presentation (2, 5, 6).

### References:

- Liu, W. *et al.* (2004) *J. Biol. Chem.* **279**:18748.
- Caminschi, I. *et al.* (2001) *Mol. Immunol.* **38**:365.
- Curtis, B.M. *et al.* (1992) *Proc. Natl. Acad. Sci. USA* **89**:8356.
- Anthony, R.M. *et al.* (2008) *Proc. Natl. Acad. Sci. USA* **105**:19571.
- Geijtenbeek, T.B. *et al.* (2000) *Cell* **100**:575.
- Garcia-Vallejo, J.J. and Y. van Kooyk (2013) *Trends Immunol.* **34**:482.