

DESCRIPTION

Species Reactivity	Mouse
Specificity	Detects mouse IGSF9B in direct ELISAs.
Source	Monoclonal Mouse IgG _{2B} Clone # 993107
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Mouse myeloma cell line NS0-derived recombinant mouse IGSF9B Arg21-Leu730 Accession # Q05BQ1
Conjugate	Alexa Fluor 647 Excitation Wavelength: 650 nm Emission Wavelength: 668 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Flow Cytometry	0.25-1 µg/10 ⁶ cells	HEK293 Human Cell Line Transfected with Mouse IGSF9 and eGFP

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. <ul style="list-style-type: none"> 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

IGSF9, also known as Dasm1, is an approximately 130 kDa transmembrane protein that plays a role in neuronal synapse maintenance and function (1). Mature mouse Dasm1 consists of a 714 amino acid (aa) extracellular domain (ECD) with 5 Ig-like domains and 2 fibronectin type-3 domains, a 21 aa transmembrane segment, and a 424 aa cytoplasmic domain (2). Within the ECD, mouse Dasm1 shares 90% and 96% aa sequence identity with human and rat Dasm1, respectively. Alternative splicing generates an additional isoform that is truncated following the fifth Ig-like domain. Dasm1 is expressed in the dorsal root and trigeminal ganglia, forebrain, cortex, dentate gyrus, pyramidal cells, Purkinje cells, and hippocampal CA1 interneurons (2-5). It localizes to dendrites, cell bodies, and post-synaptic densities (3, 6). Dasm1 functions as a homophilic adhesion protein that supports the maintenance of inhibitory synapses as well as inhibitory neurotransmission (5, 6).

References:

- Hansen, M. and P.S. Walmod (2013) *Neurochem. Res.* **38**:1236.
- Doudney, K. *et al.* (2002) *Genomics* **79**:663.
- Shi, S.-H. *et al.* (2004) *Proc. Natl. Acad. Sci. USA* **101**:13341.
- Mishra, A. *et al.* (2008) *Mol. Cell. Biol.* **28**:2782.
- Mishra, A. *et al.* (2014) *J. Neurosci.* **34**:4187.
- Shi, S.-H. *et al.* (2004) *Proc. Natl. Acad. Sci. USA* **101**:13346.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc. and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.