

## DESCRIPTION

<b>Species Reactivity</b>	Human
<b>Specificity</b>	Detects human BAI2 in direct ELISAs.
<b>Source</b>	Recombinant Monoclonal Rabbit IgG Clone # 2528B
<b>Purification</b>	Protein A or G purified from hybridoma culture supernatant
<b>Immunogen</b>	Human embryonic kidney cell, HEK293-derived human BAI2 Phe21-Asp911 Accession # O60241-3
<b>Conjugate</b>	Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm
<b>Formulation</b>	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide.  *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

## APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

	Recommended Concentration	Sample
<b>Flow Cytometry</b>	0.25-1 µg/10 <sup>6</sup> cells	HEK293 Human Cell Line Transfected with Human BAI2 and eGFP

## PREPARATION AND STORAGE

<b>Shipping</b>	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
<b>Stability &amp; Storage</b>	<b>Protect from light. Do not freeze.</b> <ul style="list-style-type: none"> <li>12 months from date of receipt, 2 to 8 °C as supplied.</li> </ul>

## BACKGROUND

Brain angiogenesis inhibitor 2 (BAI2), along with BAI1 and BAI3, form a subfamily of brain specific adhesion GPCRs. They are large, glycosylated molecules, approximately 180-220 kDa, and play a role in neuronal synapse formation and maintenance. Mature human BAI2 consists of a 903 amino acid (aa) N-terminal extracellular domain (ECD) containing four TSP1 domains, a hormone binding domain (HBD), a GAIN domain and a GPS protease sensitive linker, followed by a region with seven transmembrane segments and a 410 aa C-terminal cytoplasmic domain (1). Within the N-terminal ECD, human BAI2 shares 96% aa sequence identity with both mouse and rat BAI2. BAI2 is expressed primarily in neurons and astrocytes of the hippocampus and cerebral cortex in the brain (2, 3). To become functional, BAI2 might require activation by cleavage at the GPS domain (4). BAI family members are known to bind the C1q-like complement (C1qL) family via thrombospondin type 1 repeats (TSRs) which are conserved between BAI members (5). BAI2 has been shown to suppress expression of VEGF leading to increased neurogenesis in the dentate gyrus of the hippocampus (6).

### References:

- Duman, J.G. *et al.* (2016) *Neural Plast.* **2016**:8301737.
- Shiratsuchi, T. *et al.* (1997) *Cytogenet. Cell Genet.* **79**:103.
- Kee, H.J. *et al.* (2002) *J. Cereb. Blood Flow Metab.* **22**:1054.
- Daisuke, O. *et al.* (2010) *J. Recept Signal Transduct Res* **30**:143.
- Bolliger, M.F. *et al.* (2011) *Proc. Natl. Acad. Sci. USA* **108**:2534.
- Jeong, B.C. *et al.* (2016) *FEBS Lett* **580**:669.

## PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.