

***Affinity Purified
Goat Anti-human IgG Fc Antibody***

Catalog Number: G-102-C
Lot Numbers: WBT1519101
WBT1620081

Product Specifications

Amount	◆ 0.1 mg
Form	◆ Lyophilized
Purity	◆ Affinity purified ◆ ≥ 95% by SDS-PAGE
Buffer	◆ 10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS) with 5% trehalose.
Preservative	◆ 0.05% (w/v) sodium azide
Reconstitution	◆ Reconstitute with sterile PBS. If 0.5 mL of PBS is used, the antibody concentration will be 0.2 mg/mL.
Storage	◆ 2° - 8° C
Specificity	◆ This antibody reacts with heavy chains on human IgG, but not with light chains on immunoglobulins. No reactivity is observed to non-immunoglobulin human serum proteins by IEP.
Note	◆ Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal.