

*Affinity Purified Goat
Anti-human Creatine Kinase, Muscle (CKMM)
Antibody*

Catalog Number: G-113-C
Lot Numbers: WCK2020061
WCK2220111

Product Specifications

Amount	◆ 100 µg
Formulation	◆ Lyophilized
Purity	◆ Affinity purified
Buffer	◆ 10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 with 5% trehalose.
Preservative	◆ 0.05% (w/v) sodium azide
Reconstitution	◆ Reconstitute with sterile PBS. If 0.5 mL of PBS is used, the antibody concentration will be 0.2 mg/mL.
Storage	◆ 2° - 8° C
Note	◆ Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal.