

Affinity Purified
Goat Anti-mouse IgG Fc Antibody

Catalog Number: G-202-C
Lot Number: WDW1420111

Product Specifications

Amount	◆ 100 µg
Form	◆ Lyophilized
Purity	◆ Affinity purified
Buffer	◆ 10 mM sodium phosphate, 0.85% (w/v) NaCl, pH 7.2 (PBS) with 5% trehalose
Preservative	◆ 0.05% (w/v) sodium azide
Reconstitution	◆ Reconstitute with sterile PBS. If 0.5 mL of PBS is used, the antibody concentration will be 0.2 mg/mL.
Storage	◆ 2° - 8° C
Specificity	◆ This antibody reacts with heavy chains on mouse IgG, but not with light chains on mouse immunoglobulins. No reactivity is observed to non-immunoglobulin mouse serum proteins.
Note	◆ Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal.