

DESCRIPTION


| | |
|---------------------------|---|
| Species Reactivity | Goat |
| Specificity | Detects goat IgG heavy and light chains in direct ELISAs and Western blots. In Western blots, less than 5% cross-reactivity with mouse IgG, rabbit IgG and human IgG is observed. |
| Source | Polyclonal Chicken IgY |
| Purification | Antigen Affinity-purified from egg yolks |
| Immunogen | Goat IgG |
| Formulation | Supplied as a 0.2 µm filtered solution in Sodium Phosphate and Proclin 300. See Certificate of Analysis for details. |

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

Western Blot Use at a 1:1000 dilution.

DATA



PREPARATION AND STORAGE

Shipping The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage **Do not freeze.**

- 6 months from date of receipt, 2 to 8 °C as supplied.