

DESCRIPTION


Specificity	Detects N-terminal and C-terminal HA-tagged proteins in Western blot.
Source	Recombinant Monoclonal Rabbit IgG Clone # 1049F
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	HA Tag Peptide Accession # ABB51961
Formulation	Supplied as a solution in PBS containing a preservative. See Certificate of Analysis for details.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Western Blot	1:1000-1:3000 dilution	See Below

DATA



PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	<p>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</p> <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C, as supplied. ● 1 month, 2 to 8 °C under sterile conditions after opening. ● 6 months, -20 to -70 °C under sterile conditions after opening.