

DESCRIPTION


Species Reactivity	Human
Specificity	Detects human, mouse and rat β -Catenin in Western blots.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human β -Catenin Ala2-Leu781 Accession # P35222
Conjugate	Alexa Fluor 647 Excitation Wavelength: 650 nm Emission Wavelength: 668 nm
Formulation	Supplied in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Intracellular Staining by Flow Cytometry	5 μ L/10 ⁶ cells	See Below

DATA



PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. <ul style="list-style-type: none"> 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

β -Catenin is a 94-96 kDa member of the β -Catenin family, armadillo-repeat superfamily of intracellular molecules. Human β -Catenin is 781 amino acids (aa) in length and contains three basic regions; there is a "flexible" N-terminus (aa 1-110) that undergoes phosphorylation, O-linked glycosylation, acetylation and ubiquitination, a central region (aa 14-670) that contains 12 armadillo-like repeats and engages in protein-protein interaction, and a "flexible" C-terminus that also undergoes phosphorylation. β -Catenin is widely, if not ubiquitously expressed, being found in endothelial cells, NKT cells, cortical neurons, neural stem cells, osteoblasts, fibroblasts, skeletal muscle cells, and Tregs. β -Catenin is perhaps best known as a facilitator for Wnt signaling. In this context, and in the absence of Wnt signaling, non-nuclear or E-Cadherin-associated β -Catenin is continually synthesized and degraded via a phosphorylation-dependent program. Upon Wnt engagement of its LPR:Fzd receptor complex, phosphorylation stops, allowing for an increase in its half-life and translocation into the nucleus. Here, it binds to TCF, generating a gene-activating transcriptional complex. Notably, Wnt signaling has been described as being β -Catenin and Ca⁺⁺/non- β -Catenin based. It now appears they work in tandem, with Ca⁺⁺ mobilization facilitating the translocation of β -Catenin across the nuclear membrane. Human and mouse full-length β -Catenin shows 99% aa sequence identity.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc. and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.