

DESCRIPTION


Species Reactivity	Human
Specificity	Detects human PLZF in direct ELISAs and Western blots. In direct ELISAs, no cross-reactivity with recombinant human (rh) ZBTB38, rhZNF24, rhZNF143, rhZNF206, rhZNF281, or rhZNF423 is observed.
Source	Monoclonal Mouse IgG _{2A} Clone # 6318100
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant human PLZF Met1-Gln254 Accession # Q05516
Conjugate	PerCP (Peridinin-chlorophyll Protein Complex) Excitation Wavelength: 482 and 564 nm Emission Wavelength: 675 nm
Formulation	Supplied in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Intracellular Staining by Flow Cytometry	10 μ L/10 ⁶ cells	See Below

DATA



PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. <ul style="list-style-type: none"> ● 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

Human PLZF is a 74 kDa nuclear protein that belongs to the POK family of transcriptional repressors. It is a 673 amino acid protein that contains an N-terminal BTB domain, followed by an acidic domain, a proline-rich region and a C-terminal zinc-finger domain. PLZF forms homodimers with RAR α and LAZ3 with its zinc-finger region. Alternate splice forms exist which are tissue-specific and show a deletion of either the BTB domain, the acidic region, or the proline-rich region. Human PLZF shares 96%, 97%, 96%, and 96% amino acid identity with rat, mouse, canine, and bovine PLZF, respectively.