

Human Laminin α 1 Alexa Fluor® 594-conjugated Antibody

Monoclonal Mouse IgG_{2A} Clone # 775016

Catalog Number: IC4187T

100 μ g

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human Laminin α 1 in direct ELISAs.
Source	Monoclonal Mouse IgG _{2A} Clone # 775016
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant human Laminin α 1 Leu22-Met269 Accession # P25391
Conjugate	Alexa Fluor 594 Excitation Wavelength: 590 nm Emission Wavelength: 617 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

	Recommended Concentration	Sample
Flow Cytometry	0.25-1 μ g/10 ⁶ cells	U2OS human cell line fixed with Flow Cytometry Fixation Buffer (Catalog # FC004) and permeabilized with Flow Cytometry Permeabilization/Wash Buffer I (Catalog # FC005)

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

Laminin subunit α 1 (LAMA1) is a secreted 400 kDa extracellular matrix glycoprotein that contributes to the formation of basement membrane Laminin isoforms 1 and 3. It is one of three subunits (α , β , and γ) that interact via their coiled-coil domains to form the approximately 800 kDa cruciform, disulfide-linked, Laminin heterotrimer. The 3058 amino acid (aa) residue mature human α 1 chain contains an N-terminal Laminin VI domain (aa 18-269), followed by domains V through III containing 17 EGF-like repeats, the coiled-coil domains II and I, and five globular, Laminin G-like domains. Over aa 22-269, human Laminin α 1 shares 95% and 91% aa sequence identity with canine and mouse α 1 chain, respectively.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.