

Human GPR50 Alexa Fluor® 405-conjugated Antibody

Monoclonal Mouse IgG_{2A} Clone # 461129

Catalog Number: IC4645V
100 µg

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human GPR50. Stains human GPR50 transfectants but not irrelevant transfectants.
Source	Monoclonal Mouse IgG _{2A} Clone # 461129
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	NS0 mouse myeloma cell line transfected with human GPR50 Met1-Val617 Accession # Q13585
Conjugate	Alexa Fluor 405 Excitation Wavelength: 405 nm Emission Wavelength: 421 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

	Recommended Concentration	Sample
Intracellular Staining by Flow Cytometry	0.25-1 µg/10 ⁶ cells	A172 human glioblastoma cell line fixed with paraformaldehyde and permeabilized with saponin

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. ● 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

GPR50, also known as MTR1L, is a non-glycosylated seven-transmembrane G protein-coupled receptor that is related to the melatonin receptors MT1 and MT2. GPR50 is expressed in the hippocampus, hypothalamus, and pituitary and forms 130 kDa homodimers. It heterodimerizes with either MT1 or MT2, resulting in inhibition of MT1 but not MT2 function. An alternately spliced isoform of GPR50 has a 4 aa deletion in the large C-terminal cytoplasmic domain. The presence of this deletion as well as various polymorphisms have been associated with elevated serum triglyceride and HDL levels. The deletion may also be associated with the development of bipolar disorder. Human GPR50 shares approximately 70% amino acid sequence identity with mouse and rat GPR50.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc., and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.