

DESCRIPTION


Species Reactivity	Mouse
Specificity	Detects mouse NFIL3/E4BP4 in direct ELISAs.
Source	Recombinant Monoclonal Rabbit IgG Clone # 1218A
Purification	Protein A or G purified from cell culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant mouse NFIL3/E4BP4 Lys140-Arg462 Accession # O08750
Conjugate	Phycoerythrin Excitation Wavelength: 488 nm Emission Wavelength: 565-605 nm
Formulation	Supplied in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Intracellular Staining by Flow Cytometry	10 μ L/10 ⁶ cells	See Below

DATA



PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. <ul style="list-style-type: none"> 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

NFIL3 (Nuclear Factor, Interleukin 3 Regulated), also called E4BP4 (adenovirus E4 promoter binding protein 4), is an approximately 58 kDa transcription factor belonging to the bZIP (basic leucine zipper) family. It is a transcriptional regulator expressed in Th2 cells, Th1 cells, pro-B lymphocytes and cardiomyocytes. NFIL3 binds IL-3 and promoters as a homodimer, and enhances cell survival. NFIL3 has also been shown to be involved in the development of Innate Lymphoid Cells (ILCs). The 462 amino acid (aa) human NFIL3 contains a basic DNA binding domain (aa 73-95), a leucine zipper (aa 96-124), and a transcriptional repression domain (aa 299-363) that is potentially regulated by phosphorylation at Ser301 and Ser353. Within the region used as an immunogen, human NFIL3 shares 83% aa sequence identity with mouse and rat NFIL3.