

DESCRIPTION


Species Reactivity	Human
Specificity	Detects human Epregrulin in direct ELISAs.
Source	Monoclonal Rat IgG _{2A} Clone # 870126
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant human Epregrulin Val63-Leu108 Accession # O14944
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied either lyophilized or as a 0.2 µm filtered solution in PBS.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

ELISA	This antibody functions as an ELISA capture antibody when paired with Goat Anti-Human Epregrulin Antigen Affinity-purified Polyclonal Antibody(Catalog # AF1195). <i>This product is intended for assay development on various assay platforms requiring antibody pairs. We recommend the Human Epregrulin DuoSet ELISA Kit (Catalog # DY1195-05) for convenient development of a sandwich ELISA.</i>
--------------	---

DATA



PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.5 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> • 12 months from date of receipt, -20 to -70 °C as supplied. • 1 month, 2 to 8 °C under sterile conditions after reconstitution. • 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Epiregulin is a member of the EGF family of growth factors which includes, among others, epidermal growth factor (EGF), transforming growth factor (TGF)-alpha, amphiregulin (ARG), HB (heparin-binding)-EGF, betacellulin, and the various heregulins. All EGF family members are synthesized as transmembrane precursors and are converted to soluble forms by proteolytic cleavage. Epiregulin was originally purified from the mouse fibroblast-derived tumor cell line NIH3T3/T7 (1). The human epiregulin cDNA encodes a 169 amino acid (aa) residues transmembrane precursor with a 29 aa signal peptide, a 21 aa transmembrane domain and a 21 aa cytoplasmic domain. The putative soluble mature Epiregulin comprising the EGF-like domain (aa residues 64-104) is formed by proteolytic removal of the propeptide regions (2). There is 85% aa sequence homology between human and mouse epiregulins. Epiregulin is expressed primarily in the placenta and macrophages (3). High level expression has also been detected in various carcinomas. Epiregulin specifically binds EGFR (ErbB1) and ErbB4 but not ErbB2 and ErbB3. It activates the homodimers of both ErbB1 and ErbB4. In addition, epiregulin can also activate all possible heteromeric combinations of the four ErbB family members (4). Epiregulin stimulates the proliferation of fibroblasts, smooth muscle cells and hepatocytes. It has been shown to be an autocrine growth factor for epidermal keratinocytes as well as mesangial cells (5, 6). Epiregulin has also been shown to inhibit growth of several epithelial tumor cells. In addition, Epiregulin has been implicated in the implantation process during pregnancy (7).

References:

1. Toyoda, H. *et al.* (1995) *J. Biol. Chem.* **270**:7495.
2. Toyoda, H. *et al.* (1997) *Biochem. J.* **326**:69.
3. Komurasaki, T. *et al.* (1997) *Oncogene* **15**:2841.
4. Shelly, M. *et al.* (1998) *J. Biol. Chem.* **273**:10496.
5. Shirakata, Y. *et al.* (2000) *J. Biol. Chem.* **275**:5748.
6. Mishre, R. *et al.* (2002) *Am. J. Physiol. Renal. Physiol.* **283**:F1151.
7. Das, S.K. *et al.* (1997) *Dev. Biol.* **190**:178.