

DESCRIPTION

Species Reactivity	Mouse
Specificity	Detects mouse Frizzled-6 in direct ELISAs and Western blots. In direct ELISAs, no cross-reactivity with recombinant human Frizzled-5 or recombinant mouse Frizzled-1, -2, -3, -4, -5, -7, -8, or -9 is observed.
Source	Monoclonal Rat IgG _{2A} Clone # 211435
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Mouse myeloma cell line NS0-derived recombinant mouse Frizzled-6 His19-Val153 (predicted) Accession # Q61089
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied either lyophilized or as a 0.2 µm filtered solution in PBS.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the [Technical Information](#) section on our website.

	Recommended Concentration	Sample
Western Blot	1 µg/mL	Recombinant Mouse Frizzled-6

PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.5 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Frizzled-6 is a member of the Frizzled/Smoothed family of G-protein coupled receptors that mediate the Wnt beta-Catenin signaling pathway.