

**DESCRIPTION**


<b>Species Reactivity</b>	Human
<b>Specificity</b>	Detects recombinant human DLL1 in direct ELISAs and Western blots. In direct ELISAs and Western blots, no cross-reactivity with recombinant mouse DLL4 is observed.
<b>Source</b>	Monoclonal Mouse IgG <sub>2B</sub> Clone # 251127
<b>Purification</b>	Protein A or G purified from hybridoma culture supernatant
<b>Immunogen</b>	Mouse myeloma cell line NS0-derived recombinant human DLL1 Ser22-Gly540 Accession # AAG09716
<b>Formulation</b>	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied either lyophilized or as a 0.2 µm filtered solution in PBS.

**APPLICATIONS**

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	<b>Recommended Concentration</b>	<b>Sample</b>
<b>Western Blot</b>	1 µg/mL	Recombinant Human DLL1 (Catalog # 1818-DL) under non-reducing conditions only
<b>Flow Cytometry</b>	2.5 µg/10 <sup>6</sup> cells	See Below
<b>Human DLL1 Sandwich Immunoassay</b>		<b>Reagent</b>
<b>ELISA Capture</b>	2-8 µg/mL	Human DLL1 Antibody (Catalog # MAB1818)
<b>ELISA Detection</b>	0.1-0.4 µg/mL	Human DLL1 Biotinylated Antibody (Catalog # BAF1818)
<b>Standard</b>		Recombinant Human DLL1 (Catalog # 1818-DL)
<b>CyTOF-ready</b>	Ready to be labeled using established conjugation methods. No BSA or other carrier proteins that could interfere with conjugation.	

**DATA**



**PREPARATION AND STORAGE**

<b>Reconstitution</b>	Reconstitute at 0.5 mg/mL in sterile PBS.
<b>Shipping</b>	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
<b>Stability &amp; Storage</b>	<b>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</b> <ul style="list-style-type: none"> <li>• 12 months from date of receipt, -20 to -70 °C as supplied.</li> <li>• 1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> <li>• 6 months, -20 to -70 °C under sterile conditions after reconstitution.</li> </ul>

**BACKGROUND**

Delta-like protein 1 (DLL1) is a 90-100 kDa type I transmembrane protein that belongs to the Delta/Serrate/Lag-2 (DSL) family of Notch ligands. Mature human DLL1 consists of a 528 amino acid (aa) extracellular domain (ECD) with one DSL domain and eight EGF-like repeats, a 23 aa transmembrane segment, and a 155 aa cytoplasmic domain (1). Within the ECD, human DLL1 shares 91% aa sequence identity with mouse and rat DLL1. It shares 26%, 37%, and 54% aa sequence identity with DLL2, 3, and 4, respectively. A 60 kDa ECD fragment released by ADAM9, 12, or 17 mediated proteolysis, promotes the proliferation of hematopoietic progenitor cells (2, 3). The residual membrane-bound portion of DLL1 can be cleaved by presenilin-dependent  $\gamma$ -secretase, enabling the cytoplasmic domain to migrate to the nucleus (4). DLL1 localizes to adherens junctions on neuronal processes through its association with the scaffolding protein MAGI1 (5). DLL1 is widely expressed, and it plays an important role in embryonic somite formation, cochlear hair cell differentiation, plus B and T lymphocyte differentiation (6-11). The upregulation of DLL1 in arterial endothelial cells following injury or angiogenic stimulation is central to postnatal arteriogenesis (12). DLL1 is also overexpressed in cervical carcinoma and glioma and contributes to tumor progression (1, 13).

**References:**

1. Gray, G.E. *et al.* (1999) *Am. J. Pathol.* **154**:785.
2. Dyczynska, E. *et al.* (2007) *J. Biol. Chem.* **282**:436.
3. Karanu, F.N. *et al.* (2001) *Blood* **97**:1960.
4. Ikeuchi, T. and S.S. Sisodia (2003) *J. Biol. Chem.* **278**:7751.
5. Mizuhara, E. *et al.* (2005) *J. Biol. Chem.* **280**:26499.
6. Takahashi, Y. *et al.* (2003) *Development* **130**:4259.
7. Teppner, I. *et al.* (2007) *BMC Dev. Biol.* **7**:68.
8. Kiernan, A.E. *et al.* (2005) *Development* **132**:4353.
9. Schmitt, T.M. and J.C. Zuniga-Pflucker (2002) *Immunity* **17**:749.
10. Hozumi, K. *et al.* (2004) *Nat. Immunol.* **5**:638.
11. Santos, M.A. *et al.* (2007) *Proc. Natl. Acad. Sci.* **104**:15454.
12. Limbourg, A. *et al.* (2007) *Circ. Res.* **100**:363.
13. Purow, B.W. *et al.* (2005) *Cancer Res.* **65**:2353.