

**DESCRIPTION**


<b>Species Reactivity</b>	Human
<b>Specificity</b>	Detects human GSK-3 $\beta$ in direct ELISAs and human GSK-3 $\alpha$ and GSK-3 $\beta$ in Western blots.
<b>Source</b>	Monoclonal Mouse IgG <sub>2B</sub> Clone # 833916
<b>Purification</b>	Protein A or G purified from hybridoma culture supernatant
<b>Immunogen</b>	<i>E. coli</i> -derived recombinant human GSK-3 $\beta$ Asp200-Phe305 Accession # P49841
<b>Formulation</b>	Lyophilized from a 0.2 $\mu$ m filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied as a 0.2 $\mu$ m filtered solution in PBS.

**APPLICATIONS**

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	<b>Recommended Concentration</b>	<b>Sample</b>
<b>Western Blot</b>	2 $\mu$ g/mL	See Below

**DATA**



**PREPARATION AND STORAGE**

<b>Reconstitution</b>	Reconstitute at 0.5 mg/mL in sterile PBS.
<b>Shipping</b>	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
<b>Stability &amp; Storage</b>	<b>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</b> <ul style="list-style-type: none"> <li>● 12 months from date of receipt, -20 to -70 °C as supplied.</li> <li>● 1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> <li>● 6 months, -20 to -70 °C under sterile conditions after reconstitution.</li> </ul>

**BACKGROUND**

Glycogen Synthase Kinase-3 (GSK-3) is a serine/threonine kinase initially identified as an inhibitor of glycogen synthase. Two isoforms (GSK-3 $\alpha$  and GSK-3 $\beta$ ) share 85% amino acid identity. GSK-3 $\beta$ , inhibited by phosphorylation at S9 by Akt, is involved in energy metabolism, body pattern formation, and neuronal cell development.