

**DESCRIPTION**


<b>Species Reactivity</b>	Mouse
<b>Specificity</b>	Detects mouse ASGPR1 in direct ELISAs.
<b>Source</b>	Monoclonal Rat IgG <sub>2A</sub> Clone # 352802.113
<b>Purification</b>	Protein A or G purified from hybridoma culture supernatant
<b>Immunogen</b>	Mouse myeloma cell line NS0-derived recombinant mouse ASGPR1 Ser60-Asn284 Accession # NP_033844
<b>Formulation</b>	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied either lyophilized or as a 0.2 µm filtered solution in PBS.

**APPLICATIONS**

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

**ELISA** This antibody functions as an ELISA capture antibody when paired with Rat Anti-Mouse ASGR1/ASGPR1 Monoclonal Antibody (Catalog # [MAB27552](#)).  
*This product is intended for assay development on various assay platforms requiring antibody pairs.*

**DATA**



**PREPARATION AND STORAGE**

<b>Reconstitution</b>	Reconstitute at 0.5 mg/mL in sterile PBS.
<b>Shipping</b>	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
<b>Stability &amp; Storage</b>	<b>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</b> <ul style="list-style-type: none"> <li>● 12 months from date of receipt, -20 to -70 °C as supplied.</li> <li>● 1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> <li>● 6 months, -20 to -70 °C under sterile conditions after reconstitution.</li> </ul>

## BACKGROUND

The mouse asialoglycoprotein receptor (ASGP-R) is an endocytic recycling receptor that belongs to the long-form subfamily of the C-type/Ca<sup>++</sup>-dependent lectin family (1-3). It is a complex of two non-covalently linked subunits, a major 42 kDa glycoprotein (ASGPR1), and a minor 51 kDa glycoprotein (ASGR2). The major mouse ASGP-R subunit, ASGPR1, is synthesized as a 284 amino acid (aa) type II transmembrane (TM) protein that contains a 39 aa cytoplasmic region, a 21 aa TM segment, and a 224 aa extracellular domain (ECD) (4-6). The ECD contains two important structural regions. The first is a stalk region of 56 aa (aa's 59-117) that contributes to non-covalent oligomerization. The second is a 118 aa, carbohydrate-binding, Ca<sup>++</sup>-dependent C-type lectin domain (aa's 160-277) that is unusually stabilized by three Ca<sup>++</sup> ions (3, 5). There are two potential alternate splice forms for ASGPR1. Both are TM and show a deletion of the C-type lectin domain. One is 113 aa in length and shows a deletion of aa's 114-284 (7). The second is 132 aa in length and shows a deletion of aa's 118-146 and aa's 162-284 (8). Mouse ASGPR1 ECD is 89% and 79% aa identical to the ASGPR1 ECD in rat and human, respectively. The minor mouse ASGP-R subunit, ASGR2, is also a C-type lectin that shares the same structural organization as ASGR-1. It is 301 aa in length and has two 45 kDa and 51 kDa differentially-glycosylated isoforms (4, 6, 9). The ECD of ASGR2 is 50% aa identical to the ECD of ASGPR1. Although ASGPR1 and 2 can be expressed individually, a fully functional and stable ASGP-R requires simultaneous expression of both subunits (10-12). The stoichiometry of a functional ASGP-R is suggested to be either a 2:2, 3:1 or 3:2 ratio of ASGPR1:ASGR2 (13, 14). ASGPR1 is reported to bind Gal (nonreducing), GalNAc, and sialic acid<sub>2,6</sub>GalNAc (3, 15, 16). This is generally in the context of triantennary or tetraantennary configurations (2).

## References:

1. Stockert, R. J. (1995) *Physiol. Rev.* **75**:591.
2. Weigel, P.H. and J.H.N. Yik (2002) *Biochim. Biophys. Acta* **1572**:341.
3. Meier, M. *et al.* (2000) *J. Mol. Biol.* **300**:857.
4. Takezawa, R. *et al.* (1993) *Biochim. Biophys. Acta* **1172** :220.
5. Monroe, R.S. and B.E. Huber (1994) *Gene* **148**:237.
6. Sanford, J.P. *et al.* (1988) *DNA* **7**:721.
7. Heath, P. (2005) GenBank Accession #:Q5NCV2
8. Heath, P. (2005) GenBank Accession #:Q5NCV1.
9. Sanford, J.P. and D. Doyle (1990) *Biochim. Biophys. Acta* **1087**:259.
10. Braun, J.R. *et al.* (1996) *J. Biol. Chem.* **271**:21160.
11. Ishibashi, S. *et al.* (1994) *J. Biol. Chem.* **269**:27803.
12. Tozawa, R. *et al.* (2001) *J. Biol. Chem.* **276**:12624.
13. Bider, M.D. *et al.* (1996) *J. Biol. Chem.* **271**:31996.
14. Lodish, H. (1991) *Trends Biochem. Sci.* **16**:374.
15. Westerlind, U. *et al.* (2004) *Glyconj. J.* **21**:227.
16. Park, E.I. *et al.* (2005) *Proc. Natl. Acad. Sci. USA* **102**:17125.