

**DESCRIPTION**

<b>Species Reactivity</b>	Human
<b>Specificity</b>	Detects human NELL1 in direct ELISAs. In direct ELISAs, no cross-reactivity with recombinant human NELL2 or recombinant mouse NELL1 is observed.
<b>Source</b>	Monoclonal Mouse IgG <sub>2B</sub> Clone # 756435
<b>Purification</b>	Protein A or G purified from hybridoma culture supernatant
<b>Immunogen</b>	Chinese hamster ovary cell line CHO-derived recombinant human NELL1 Arg17-Asn810 Accession # Q92832
<b>Formulation</b>	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied as a 0.2 µm filtered solution in PBS.


**APPLICATIONS**

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	<b>Recommended Concentration</b>	<b>Sample</b>
<b>Immunocytochemistry</b>	8-25 µg/mL	See Below

**DATA**

**Immunocytochemistry**



**NELL1 in IMR32 Human Cell Line.** NELL1 was detected in immersion fixed IMR32 human neuroblastoma cell line using Mouse Anti-Human NELL1 Monoclonal Antibody (Catalog # MAB54871) at 10 µg/mL for 3 hours at room temperature. Cells were stained using the NorthernLights™ 557-conjugated Anti-Mouse IgG Secondary Antibody (red; Catalog # NL007) and counterstained with DAPI (blue). Specific staining was localized to cytoplasm. View our protocol for [Fluorescent ICC Staining of Cells on Coverslips](#).

**PREPARATION AND STORAGE**

<b>Reconstitution</b>	Sterile PBS to a final concentration of 0.5 mg/mL.
<b>Shipping</b>	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
<b>Stability &amp; Storage</b>	<b>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</b> <ul style="list-style-type: none"> <li>● 12 months from date of receipt, -20 to -70 °C as supplied.</li> <li>● 1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> <li>● 6 months, -20 to -70 °C under sterile conditions after reconstitution.</li> </ul>

**BACKGROUND**

NELL1 (neural EGF-like like protein 1) is an 89 kDa (predicted) member of the EGF-like domain containing family, Laminin G/N-TSP1/Pentraxin gene superfamily of molecules. When secreted, NELL1 exists as a phosphoglycoprotein that can add as much as 50 kDa to the calculated MW. NELL1 has restricted expression, being limited to pre-B cells and osteoblasts, where it apparently promotes osteoblast maturation and bone formation. In tumors, it is found in neuroblastoma-derived cells. NELL1 is both secreted and retained intracellularly where it is phosphorylated by PKC. The human NELL1 precursor is 810 amino acids (aa) in length. It contains a 16 aa signal sequence plus a 794 aa mature region. The mature region possesses an N-terminal TSP domain (aa 81-230), two VWFC domains (aa 271-390), six consecutive EGF-like domains (aa 391-631), and three additional C-terminal VWFC domains (aa 632-807). Secreted NELL1 forms a 400-420 kDa noncovalent homotrimer. Over aa 17-810, human NELL1 shares 93% aa identity with mouse and rat NELL1. Alternate splicing generates an additional isoform of human NELL1 that lacks the fifth EGF-like domain.