

DESCRIPTION


Species Reactivity	Human
Specificity	Detects human Sphingosine Kinase 1/SPHK1 in direct ELISAs.
Source	Monoclonal Mouse IgG ₁ Clone # 647706
Purification	Protein A or G purified from cell culture supernatant
Immunogen	<i>S. frugiperda</i> insect ovarian cell line Sf 21-derived recombinant human Sphingosine Kinase 1/SPHK1 Asp2-Leu398 Accession # Q9NYA1
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied either lyophilized or as a 0.2 µm filtered solution in PBS.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

ELISA	This antibody functions as an ELISA capture antibody when paired with Sheep Anti-Human Sphingosine Kinase 1/SPHK1 Antigen Affinity-purified Polyclonal Antibody (Catalog # AF5536). <i>This product is intended for assay development on various assay platforms requiring antibody pairs. We recommend the Human Sphingosine Kinase 1/SPHK1 DuoSet ELISA Kit (Catalog # DY5536) for convenient development of a sandwich ELISA.</i>
--------------	--

DATA



PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.5 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Sphingosine kinases are cytosolic or membrane-associated enzymes that catalyze the phosphorylation of sphingosine to sphingosine-1-phosphate (S1P). Two types of sphingosine kinases, SPHK1 and SPHK2, are known to be expressed in human cells. The two enzymes share considerable amino acid sequence similarity, but differ in their N-terminal and central regions (1). The two proteins also differ in tissue distribution and some kinetic properties (1). S1P is a lipid messenger that regulates diverse physiological processes including cell proliferation, migration, apoptosis, inflammation, calcium homeostasis and cytoskeletal structure (2, 3). The level of S1P is tightly controlled by SPHKs and S1P degrading enzymes. SPHK1 and its activation can be stimulated by several growth factors such as tumor necrosis factor- α , epidermal growth factor and transforming growth factor- β (3, 4). Expression of SPHK1 has been found to increase in many human solid tumors and overexpression of SPHK1 is associated with tumor angiogenesis (5). Such studies have implicated SPHK1 as a new target for cancer treatment.

References:

1. Liu, H. *et al.* (2000) *J. Biol. Chem.* **275**:19513.
2. Spiegel, S. (1999) *J. Leukocyte Biol.* **65**:341.
3. Alemany, R. *et al.* (2007) *Naunyn-Schmiedeberg's Arch. Pharmacol.* **374**:413.
4. Pederson, L. *et al.* (2008) *Proc. Natl. Acad. Scis USA.* **105**:20764.
5. Shida, D. *et al.* (2008) *Curr. Drug Targets.* **9**:662.