

## DESCRIPTION



<b>Species Reactivity</b>	Human
<b>Specificity</b>	Detects human Proprotein Convertase 2/PCSK2 in direct ELISAs. In direct ELISAs and Western blots, no cross-reactivity with recombinant human (rh) Proprotein Convertase 1, 7, 9, or rhFurin is observed.
<b>Source</b>	Monoclonal Mouse IgG <sub>2B</sub> Clone # 694009
<b>Purification</b>	Protein A or G purified from hybridoma culture supernatant
<b>Immunogen</b>	Chinese hamster ovary cell line CHO-derived recombinant human Proprotein Convertase 2/PCSK2 Gly110-Asn638 Accession # P16519
<b>Formulation</b>	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied as a 0.2 µm filtered solution in PBS.

## APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
<b>Western Blot</b>	2 µg/mL	See Below
<b>Immunocytochemistry</b>	8-25 µg/mL	See Below
<b>Immunoprecipitation</b>	25 µg/mL	Conditioned cell culture medium spiked with Recombinant Human Proprotein Convertase 2/PCSK2 (Catalog # 6018-SE), see our available <a href="#">Western blot detection antibodies</a>

## DATA

<p><b>Western Blot</b></p>  <p><b>Detection of Human Proprotein Convertase 2/PCSK2 by Western Blot.</b> Western blot shows lysates of human brain (cerebellum) tissue. PVDF Membrane was probed with 2 µg/mL of Human Proprotein Convertase 2/PCSK2 Monoclonal Antibody (Catalog # MAB6018) followed by HRP-conjugated Anti-Mouse IgG Secondary Antibody (Catalog # HAF007). A specific band was detected for Proprotein Convertase 2/PCSK2 at approximately 75 kDa (as indicated). This experiment was conducted under reducing conditions and using Immunoblot Buffer Group 1.</p>	<p><b>Immunocytochemistry</b></p>  <p><b>Protein Convertase 2/PCSK2 in HepG2 Human Cell Line.</b> Proprotein Convertase 2/PCSK2 was detected in immersion fixed HepG2 human hepatocellular carcinoma cell line using Mouse Anti-Human Proprotein Convertase 2/PCSK2 Monoclonal Antibody (Catalog # MAB6018) at 10 µg/mL for 3 hours at room temperature. Cells were stained using the NorthernLights™ 557-conjugated Anti-Mouse IgG Secondary Antibody (red; Catalog # NL007) and counterstained with DAPI (blue). View our protocol for <a href="#">Fluorescent ICC Staining of Cells on Coverslips</a>.</p>
--	--

## PREPARATION AND STORAGE

<b>Reconstitution</b>	Sterile PBS to a final concentration of 0.5 mg/mL.
<b>Shipping</b>	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
<b>Stability &amp; Storage</b>	<b>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</b> <ul style="list-style-type: none"> <li>● 12 months from date of receipt, -20 to -70 °C as supplied.</li> <li>● 1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> <li>● 6 months, -20 to -70 °C under sterile conditions after reconstitution.</li> </ul>

## BACKGROUND

Proprotein Convertase 2 (PCSK2) is a subtilisin-like serine peptidase that processes proteins into biologically active products in the secretory pathway. Like other furin/kexin family members, PCSK2 cleaves its substrates primarily after sites of paired basic amino acid residues (1). Autoactivation of PCSK2 occurs in a post-Golgi compartment of the secretory system (2). Interaction with the secretory protein 7B2 is required for proper activation of PCSK2 (3). The N-terminal domain of 7B2 stabilizes active PCSK2, and a C-terminal fragment of 7B2 can inhibit PCSK2 activity (4). rhPCSK2 was coexpressed with 7B2 to facilitate activation of the enzyme. PCSK2 is a major proteolytic processing enzyme of the regulated secretory pathway of the neuroendocrine system, where it generates a number of hormones and neuropeptides (5). Products of PCSK2 processing include enkephalins, insulin, somatostatin, dynorphin, and LHRH. Within the region used as an immunogen, human PCSK2 shares 97% aa identity with mouse PCSK2.

### References:

1. Zhou, A. *et al.* (1999) *J. Biol. Chem.* **274**:20745.
2. Lamango, N.S. *et al.* (1999) *Arch. Biochem. Biophys.* **362**:275.
3. Zhu, X. and I. Lindberg (1995) *J. Cell Biol.* **129**:1641.
4. Lamango, N. *et al.* (1996) *Arch. Biochem. Biophys.* **330**:238.
5. Rouillé, Y. *et al.* (1995) *Front. Neuroendocrinol.* **16**:322.