

DESCRIPTION


Species Reactivity	Human
Specificity	Detects human TRIM21 in direct ELISAs and Western blots. In direct ELISAs, no cross-reactivity with recombinant human (rh) TRIM5 or rhTRIM32 is observed.
Source	Monoclonal Mouse IgG _{2A} Clone # 672714
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant human TRIM21 Arg195-Pro293 Accession # P19474
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied as a 0.2 µm filtered solution in PBS.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Western Blot	1 µg/mL	See Below

DATA



PREPARATION AND STORAGE

Reconstitution	Sterile PBS to a final concentration of 0.5 mg/mL.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

TRIM21 (Tripartite motif-containing protein 21; also Ro(SS-A), 52 kDa Ro Protein/Ro52, and RING finger protein 81) is a 52-56 kDa member of the RING finger-B-box-coiled-coil family of proteins. It is an E3 ligase that is found in both nucleus and cytoplasm, where it is often associated with microtubules. TRIM21 ubiquitinates select proteins. In B cells, it targets the Fc fragment of misfolded IgG, providing QC on its production. In macrophages, it acts in a nondegradative manner on IRF8, promoting innate immunity. Human TRIM21 is 475 amino acids (aa) in length. It contains one E3 ligase RING finger domain (aa 16-55), a B-Box type zinc-finger region (aa 92-123), a coiled-coil region (aa 128-238) and a C-terminal SPRY/B30.2 Ig binding domain (aa 268-465). TRIM21 is reported to form trimers. Over aa 195-293, human TRIM21 exhibits 72% aa identity with mouse TRIM21.