

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human LIGHT/TNFSF14 in direct ELISAs.
Source	Recombinant Monoclonal Mouse IgG ₁ Clone # 115515R
Purification	Protein A or G purified from cell culture supernatant
Immunogen	Mouse myeloma cell line NS0-derived recombinant human LIGHT/TNFSF14 Asp74-Val240 Accession # O46557
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied either lyophilized or as a 0.2 µm filtered solution in PBS.


APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

ELISA This antibody functions as an ELISA capture antibody when paired with Goat Anti-Human LIGHT/TNFSF14 Antigen Affinity-purified Polyclonal Antibody (Catalog # [AF664](#)).

This product is intended for assay development on various assay platforms requiring antibody pairs. We recommend the Human LIGHT/TNFSF14 DuoSet ELISA Kit (Catalog # [DY664](#)) for convenient development of a sandwich ELISA or the Human LIGHT/TNFSF14 Quantikine ELISA Kit (Catalog # [DLIT00](#)) for a complete optimized ELISA.

DATA



PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.5 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Human LIGHT, also known as TNFSF14, is a type II membrane protein that is a member of the TNF superfamily. LIGHT is an acronym which stands for "is homologous to lymphotoxins, exhibits inducible expression, and competes with HSV glycoprotein D for HVEM, a receptor expressed by T lymphocytes". LIGHT has also been called HVEM-L and LT- γ . LIGHT is a 240 amino acid (aa) protein that contains a 37 aa cytoplasmic domain, a 22 aa transmembrane region, and a 181 aa extracellular domain. Similar to other TNF ligand family members, LIGHT is predicted to assemble as a homotrimer. LIGHT is produced by activated T cells and was first identified by its ability to compete with HSV glycoprotein D for HVEM binding. LIGHT has also been shown to bind to the lymphotoxin beta receptor (LT β R) and the decoy receptor (DcR3/TR6). LIGHT overexpression in tumor cells induces apoptosis, which can be enhanced by IFN- γ .

References:

1. Mauri, D.N. *et al.* (1998) *Immunity* **8**:21.
2. Zhai, Y. *et al.* (1998) *J. Clin. Invest.* **102**:1142.
3. Harrop, J.A. *et al.* (1998) *J. Biol. Chem.* **273**:27548.
4. Yu, K-Y. *et al.* (1999) *J. Biol. Chem.* **274**:13733.