

DESCRIPTION


Species Reactivity	Mouse
Specificity	Detects mouse CCL27/CTACK in direct ELISAs and Western blots. In Western blots, this antibody does not cross-react with recombinant cotton rat CCL5, recombinant human CCL1, 2, 3, 4, 5, 7, 8, 11, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 25, 27, 28, recombinant mouse CCL1, 2, 3, 4, 6, 7, 9/10/MIP-1 γ, 11, 12, 19, 21, 22, 25, or recombinant rat CCL20.
Source	Monoclonal Rat IgG _{2B} Clone # 68623
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant human CCL27/CTACK Leu26-Asn120 Accession # NP_035466
Endotoxin Level	<0.10 EU per 1 μg of the antibody by the LAL method.
Formulation	Lyophilized from a 0.2 μm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied as a 0.2 μm filtered solution in PBS.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

Mouse CCL27/CTACK Sandwich Immunoassay	Reagent
ELISA Capture	2-8 μg/mL
ELISA Detection	0.1-0.4 μg/mL
Standard	Mouse CCL27/CTACK Antibody (Catalog # MAB725) Mouse CCL27/CTACK Biotinylated Antibody (Catalog # BAF725) Recombinant Mouse CCL27/CTACK (Catalog # 725-CK)
Neutralization	Measured by its ability to neutralize CCL27/CTACK-induced chemotaxis in the BaF3 mouse pro-B cell line transfected with mouse CCR10. The Neutralization Dose (ND ₅₀) is typically 2.5-10 μg/mL in the presence of 1 μg/mL Recombinant Human CCL27/CTACK.

DATA



PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.5 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

CCL27, also known as CTACK (cutaneous T cell-attracting chemokine), ALP, ILC, and ESkin, is a member of the CC family of chemokines (1). Mature mouse CCL27 is a 95 amino acid (aa) protein that shares 57% aa and 87% aa sequence identity with human and rat CCL27, respectively (2). It shares 18-31% aa sequence identity with other mouse CC chemokines. An alternately spliced form of mouse CCL27, known as PESKY, is localized to the nucleus and promotes cellular migration (3). CCL27 is constitutively expressed by keratinocytes and is upregulated by inflammatory stimuli and in wounded skin (4-7). CCL27 binds the chemokine receptor CCR10, glycosaminoglycans in the extracellular matrix, sulfated tyrosine residues on PSGL-1, and determinants on the surface of fibroblasts and endothelial cells (5, 7-9). CCL27 cooperates with CCL17/TARC in inducing the migration of cutaneous lymphocyte antigen (CLA) positive memory T cells to the skin during inflammation (4, 6, 10-12). Endothelial cell-bound CCL27 can mediate the adhesion of those cells to CLA⁺ T cells (6). CCL27 also induces the migration of keratinocyte precursors from bone marrow to the skin, thereby promoting wound healing (7). In humans, serum CCL27 levels are elevated and correlate with disease severity in atopic dermatitis, psoriasis vulgaris, and mycosis fungoides (13-15).

References:

1. Kunkel, L. and E.C. Butcher (2002) *Immunity* **16**:1.
2. Hromas, R. *et al.* (1999) *Biochem. Biophys. Res. Commun.* **258**:737.
3. Gortz, A. *et al.* (2002) *J. Immunol.* **169**:1387.
4. Morales, J. *et al.* (1999) *Proc. Natl. Acad. Sci. USA* **96**:14470.
5. Homey, B. *et al.* (2000) *J. Immunol.* **164**:3465.
6. Homey, B. *et al.* (2002) *Nat. Med.* **8**:157.
7. Inokuma, D. *et al.* (2006) *Stem Cells* **24**:2810.
8. Jarmin, D. *et al.* (2000) *J. Immunol.* **164**:3460.
9. Hirata, T. *et al.* (2004) *J. Biol. Chem.* **279**:51775.
10. Vestergaard, C. *et al.* (2004) *Exp. Dermatol.* **13**:551.
11. Reiss, Y. *et al.* (2001) *J. Exp. Med.* **194**:1541.
12. Soler, D. *et al.* (2003) *Blood* **101**:1677.
13. Kakinuma, T. *et al.* (2003) *J. Allergy Clin. Immunol.* **111**:592.
14. Hijnen, D. *et al.* (2004) *J. Allergy Clin. Immunol.* **113**:334.
15. Fujita, Y. *et al.* (2006) *Clin. Cancer Res.* **12**:2670.