

DESCRIPTION


Species Reactivity	Human
Specificity	Detects human XIAP in direct ELISAs.
Source	Monoclonal Mouse IgG _{2B} Clone # 117305
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant human XIAP Met1-Ser497 Accession # P98170
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied either lyophilized or as a 0.2 µm filtered solution in PBS.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

ELISA	This antibody functions as an ELISA capture antibody when paired with Mouse Anti-Human XIAP Monoclonal Antibody (Catalog # MAB8223). This product is intended for assay development on various assay platforms requiring antibody pairs. We recommend the Human Total XIAP DuoSet IC ELISA Kit (Catalog # DYC822-2) for convenient development of a sandwich ELISA.
--------------	--

DATA



PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.5 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

XIAP (X-chromosome linked inhibitor of apoptosis) is a member of the apoptosis (IAP) family of proteins that inhibit caspases. The BIR2 domain of XIAP inhibits Caspase-3 and Caspase-7. The ability of XIAP to inhibit caspases is prevented by SMAC/Diablo through binding to XIAP-BIR2 and -BIR3 domains.