

DESCRIPTION


Species Reactivity	Human
Specificity	Detects human Troponin I/TNNI3 in direct ELISAs.
Source	Monoclonal Mouse IgG _{2B} Clone # 923945
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	<i>E. coli</i> -derived recombinant human Troponin I/TNNI3 Met1-Ser210 Accession # P19429
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied either lyophilized or as a 0.2 µm filtered solution in PBS.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

ELISA	This antibody functions as an ELISA detection antibody when paired with Mouse Anti-Human Troponin I/TNNI3 Monoclonal Antibody (Catalog # MAB85942). This product is intended for assay development on various assay platforms requiring antibody pairs.
--------------	--

DATA



PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.5 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> • 12 months from date of receipt, -20 to -70 °C as supplied. • 1 month, 2 to 8 °C under sterile conditions after reconstitution. • 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Troponin I, also known as TNI, is a 24-29 kDa component of a protein complex on striated muscle thin filaments. Troponin I inhibits the calcium-dependent muscle contraction mediated by Troponins C and T. The expression of cardiac Troponin I (TNNI3) is restricted to cardiac muscle, while TNNI1 and TNNI2 (encoded by distinct genes) are expressed in skeletal muscle. Mutations of cardiac Troponin I are associated with hereditary cardiomyopathy. Human cardiac Troponin I shares 93% amino acid sequence identity with mouse and rat cardiac Troponin I.