
MATERIAL DATA SHEET

**Recombinant Human Poly-Ubiquitin Wild-type Chains (2-7) (K63), Biotin
Cat. # UCB-330**

Poly-Ubiquitin chains are composed of Ubiquitin monomers that are covalently linked through isopeptide bonds, which typically form between a lysine residue of one Ubiquitin molecule and the C-terminal glycine residue of another Ubiquitin molecule (1). Each human Ubiquitin monomer is 76 amino acids (aa) in length and shares 96% and 100% aa identity with yeast and mouse Ubiquitin, respectively (2). Seven of the 76 aa in Ubiquitin are lysine residues that can participate in poly-Ubiquitin chain formation. Linkage through specific lysine residues is thought to serve as a signal that affects protein degradation, signaling, trafficking, and other cellular processes (3-8).

This mixture of poly-Ubiquitin chains contains di-Ubiquitin and higher MW species; mono-Ubiquitin has been removed. These chains have been modified with biotin via primary amine coupling. This results in multiple biotinylated species modified at the N-terminus, as well as lysine residues. Biotinylated Ubiquitin can be detected using avidin-linked reagents.

Product Information

Quantity:	50 µg
MW:	17 kDa (Ub2), 26 kDa (Ub3), 34 kDa (Ub4), 43 kDa (Ub5), 52 kDa (Ub6), and 60 kDa (Ub7)
Source:	<i>E. coli</i> -derived Accession # P0CG47
Stock:	X mg/ml in 50 mM HEPES pH 8.0
Solubility:	
Purity:	>95%, by SDS-PAGE under reducing conditions and visualized by Colloidal Coomassie® Blue stain.

Use & Storage

Use: Ubiquitin chains vary in length, linkage, and function. K63-linked Biotinylated Poly-Ubiquitin Chains (Ub2-7) are ideal for use in assays that utilize avidin-linked reagents for visualization or quantitation. Reaction conditions will need to be optimized for each specific application. **IMPORTANT:** Heating this product in SDS-PAGE buffer or terminating reactions containing this product with heated SDS-PAGE buffer could lead to unexpected, high apparent molecular weight banding or smearing on gels that is not representative of product purity. For optimal results, we recommend incubation in SDS-PAGE buffer + DTT at <40 °C for 20 minutes prior to gel electrophoresis.

Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

- 12 months from date of receipt, -20 to -70 °C as supplied.
- 3 months, -20 to -70 °C under sterile conditions after opening.

Literature

References:

1. Scheffner, M. *et al.* (1995) *Nature* **373**:81.
2. Sharp, P.M. & W.-H. Li (1987) *Trends Ecol. Evol.* **2**:328.
3. Behrends, C. & J.W. Harper (2011) *Nat. Struct. Mol. Biol.* **18**:520.
4. Greene, W. *et al.* (2012) *PLoS Pathog.* **8**:e1002703.
5. Henry, A.G. *et al.* (2012) *Dev. Cell* **23**:519.
6. Tong, X. *et al.* (2012) *J. Biol. Chem.* **287**:25280.
7. Wei, W. *et al.* (2004) *Nature* **428**:194.
8. Zhang, J. *et al.* (2012) *J. Biol. Chem.* **287**:28646.

For research use only. Not for use in humans.